

Author(s): YMa  
Date: June 2019


# Data Backup, Deletion and Restoration in EasyOne Air

## 1 Database Backup


Before deleting any data (described in Chapter 2), it is recommended to either synchronize the data with EasyOne Connect on your workstation, or back up the data onto an SD card, which is described here.

### Prerequisite:

An SD card is inserted into the EasyOne Air device. The SD card slot is located inside the battery compartment, which may be opened by turning the battery lock to the right with the help of a coin.



### Procedure:



### Important notes:

- Data backup can only be performed by users with Admin level rights, or the default Admin user if user management is not activated.
- The backup file is located in a backup folder on the SD card and named using the date/time of the backup: YYYY-MM-DD\_HH-MM-SS\_eon.db
- Should you choose to encrypt the backup, which is recommended, remember your password or save it in a secure location, as you will need it later to restore the data (described in Chapter 3).
- For security reasons, always remove the SD card from the device after performing a backup. Keep it in a safe place to avoid unauthorized access to the data. This is especially important when the backup is unencrypted.

## 2 Delete Data (Irreversible Action!)

With EasyOne Air firmware version 1.11.0 and above, you have the options to:

- Delete selected or all patient profiles and associated test data
- Delete individual test sessions
- Delete the entire database


### 2.1 Delete Patient Profile(s)


Deleting a patient profile removes all associated test data including demographic information from the internal storage of EasyOne Air.

#### Prerequisite:

Patient data is either synchronized with EasyOne Connect on your workstation or backed up in an SD card, as described in Chapter 1.

#### Procedure:





Confirm

### Important notes:

- **The deletion is permanent and cannot be undone.**
- Patient deletion can only be performed by users with Admin level rights, or the default Admin user if user management is not activated.

## 2.2 Delete Individual Test Session


This action is useful when you want to delete test data but keep the patient demographic profile.

### Procedure:



Select patient

Select test session




Press  to delete

Confirm

Admin user password

### Alternative Procedure:

Use this procedure if the session to be deleted was done recently, as the test reports are sorted by date (by default).



### Important notes:

- **The deletion is permanent and cannot be undone.**
- Test deletion can only be performed by users with Admin level rights, or the default Admin user if user management is not activated.


## 2.3 Delete Entire Database

Deleting the database permanently erases all patient profiles and associated test data. It also performs a factory reset, meaning that all your configuration settings are reset to the factory defaults.

### Prerequisite:

Patient data is either synchronized with EasyOne Connect on your workstation or backed up to an SD card, as described in Chapter 1.

### Procedure:



### Important notes:

- **The deletion is permanent and cannot be undone. It is strongly recommended to perform a backup before deleting the database.**
- All configuration settings are reset to the factory defaults.
- A database deletion can only be performed by users with Admin level rights, or the default Admin user if user management is not activated.


### 3 Restore Database from a Backup

Use this feature to return your EasyOne Air to a previous status, which is captured by the backup process described in Chapter 1.

**Prerequisite:**

The backup data is in the form of a .db file and is located within the backup folder on the SD card. Insert the SD card into the EasyOne Air device. If the backup data is encrypted, the encryption password is needed to perform the database restore.

**Procedure:**



**Important notes:**

- **All data or setting changes occurring during the time between backup and restore will be permanently lost. It is strongly recommended to perform a backup before restoring the database.**
- A database restore can only be performed by users with Admin level rights, or the default Admin user if user management is not activated.