

EasyOne FlowTube

Breathing Tube for
Infection Control

**Avoid cross-contamination
and reduce cleaning to an
absolute minimum with
the single-patient-use
EasyOne FlowTube.**



All parts exposed to patient's breath
are **single-use** and **disposable**

No expiration date

No special storage conditions

Technical data

Material

EasyOne FlowTube: Polyethylene (PP)
Bag: Polyethylene (PE)

Transport and storage information

Temp. -20 to 50°C
Rel. Humidity 5 to 90%
Athmosph. Pressure 500 to 1060 hPa
no expiration date

MDR (EU) 2017/745

Class I

Dimension (single EasyOne FlowTube)

21 x 38 x 150 mm (H x W x D), 8 g

Compatible with

EasyOne Air

Order information

Part Number

Description

5050-50	EasyOne FlowTube, individually wrapped, dispenser box of 50 pcs
5050-200	EasyOne FlowTube, individually wrapped, box of 200 pcs (dispenser box included)
5050-500	EasyOne FlowTube, individually wrapped, box of 500 pcs (without dispenser box)