

Recommended Printers for EasyOne Pro

1 Scope

This document lists the recommended printers for EasyOne Pro. The pre-installed printers listed in chapter 2 do not require any installation on the device. The driver installations for the printer drivers listed in chapter 3 and 4 can be downloaded from <http://www.ndd.ch/index/downloads-software>. For instruction on how to install the driver on a device, refer to the Application Note “Using Printers with EasyOne Pro” (<http://www.ndd.ch/index/downloads-application>).

2 Pre-Installed Printers on EasyOne Pro

- HP Universal Printer PCL5 (supports most of the printers connected on a network)
- HP OfficeJet H470 (on all EasyOne Pro images)
- HP DeskJet 6940 (on all EasyOne Pro images except V03.02 nnd A)
- HP OfficeJet 6100 (since EasyOne Pro image V03.02 nnd F)
- HP OfficeJet Pro 8100 (since EasyOne Pro image V03.02 nnd F)

Remark: The EasyOne Pro software image version can be determined via Utilities / Configuration / About / Details.

3 Recommended Printers (Tested)

- Epson Series N11, NX120, S22, SX125, T12, T13, T22, T22E, TX120, TX121
- Epson WorkForce WF-2530
- HP DeskJet 3000
- HP Inkjet 1000
- HP OfficeJet 6000

4 Printers (Untested; but drivers available)

Brother	Canon	Epson	Kyocera
HL-2140/2150N /2170W	i-SENSYS LBP6020	Inkjet L100, L200	FS-1030D
HL-2250DN	PIXMA iP1800	Series T10, T11	Universal Driver
HL-5340D	PIXMA iP2700	Stylus C79	
MFC-J435W	PIXMA iP4900	XP-312	
	PIXMA iP6220D		
	PIXMA iP90		
	PIXMA MX452		
	PIXMA MX850		

HP	Lanier	Lexmark	Samsung
Business InkJet 1200	LP020C	C543	ML-1210
DeskJet 1000		Universal Driver	ML-1750
DeskJet 1015			SCX-4200
DeskJet 2515			Universal Driver
DeskJet 3000			
DeskJet 460			
DeskJet 5440			
DeskJet 5940			
DeskJet 6520			
DeskJet 6980			
DeskJet 960			
DeskJet F2400			
InkJet 1000			
LaserJet 1018			
LaserJet 1020, 1022 PLUS			
LaserJet 1320			
LaserJet 1566			
LaserJet 4700			
LaserJet M1005			
LaserJet P1005			
LaserJet P1100,P1560, P1600 XPS			
LaserJet P2014			
LaserJet P2035			
LaserJet P2055			
LaserJet P3015			
LaserJet Pro 300/400 Color M351/M451			
OfficeJet 100			
OfficeJet 150			
OfficeJet 6000			
OfficeJet 6100			
OfficeJet J4580			
OfficeJet J4660			
OfficeJet Pro 8000			
OfficeJet Pro 8100			
Series D1400, 2400 ,4200			
Series D5500			