

Spirette™

Breathing Tube for Infection Control

Avoid cross-contamination and reduce cleaning to an absolute minimum with the single-patient-use Spirette.



All parts exposed to patient's breath are **single-use** and **disposable**

No expiration date

No special storage conditions

Technical data

Material	Spirette: Linear low density polyethylene (LLDPE) Bag: Polyethylene (PE)
Transport and storage information	Temp. -20 to 50°C Rel. Humidity 5 to 95% Atmosph. Pressure 500 to 1060 hPa no expiration date
MDR (EU) 2017/745 // MDD 93/42/EEC	Class IIa
Dimension (single spirette)	25 x 38 x 140 mm (H x W x D), 8 g
Compatible with	EasyOne, Easy on-PC, EasyOne Pro and EasyOne Pro LAB

Order information

Part Number	Description
2050-1	Spirette, individually wrapped, dispenser box of 50 pcs
2050-5	Spirette, individually wrapped, box of 200 pcs (dispenser box included)